SERION ELISA classic

ESR1201A

BORDETELLA PERTUSSIS TOXIN IgA

SGI.FE

Qualitätskontrollzertifikat / Quality Control Certificate

Kitcharge / Lot SGI.FE IFU-Version 1201-9

Verw. bis / Exp. 2020-07

06.08.2018

Prüfdatum /

Date of control



Verwendete Reagenzien / Reagents used	Standard		Standard Kur	Standard Kurve / Standard curve							
Teststreifen / Antigen coated strips	SFI.BR	Ref Werte / Ref	f. Values	5	Gültigke	itsbereich	/ Valid	dity Range	Parameter	Α	0,090
Standardserum / Standard serum	SGI.AG	OD	0,86		OD	0,43	-	1,46		В	1,471
Negativ Kontrolle / Negative control	SGI.AF									С	3,881
Konjugat / Conjugate	SGI.BP++	Units	19,6	IU/ml						D	3,783
Quantifizierungsgrenzen / Limits of quantifi	IU/mI	10	-	300							
Grenzwertbereich / Borderline range	IU/ml	15	-	20							

OD Bereich / OD Range 405 nm, Standardserum / Standard serum																				
0,43	3 -	0,47	0,48	- 0),53	0,54	- 0,58	0,59	- 0,64	0,65	- 0,69	0,70	- 0,74	0,75	- 0,80	0,81	- 0,85	0,86	IU/ml	Interpretation
	<	0,34		< 0),38		< 0,42		< 0,46		< 0,50		< 0,54		< 0,58		< 0,62	< 0,65	< 15,0	neg
0,34	- ا	0,46	0,38	- 0),52	0,42	- 0,57	0,46	- 0,63	0,50	- 0,68	0,54	- 0,74	0,58	- 0,79	0,62	- 0,85	0,65 - 0,88	15,0 - 20,0	gw / borderline
	>	0,46		> 0),52		> 0,57		> 0,63		> 0,68		> 0,74		> 0,79		> 0,85	> 0,88	> 20,0	pos

OD Bereich / OD Range 405 nm, Standardserum / Standard serum												
IU/ml	0,86	0,87 - 0,94	0,95 - 1,01	1,02 - 1,09	1,10 - 1,16	1,17 - 1,24	1,25 - 1,31	1,32 - 1,39 1,4	40 - 1,46	Interpretation		
< 15,0	< 0,65	< 0,68	< 0,74	< 0,79	< 0,85	< 0,91	< 0,96	< 1,02	< 1,08	neg		
15,0 - 20,0	0,65 - 0,88	0,68 - 0,92	0,74 - 1,00	0,79 - 1,08	0,85 - 1,15	0,91 - 1,23	0,96 - 1,31	1,02 - 1,38 1,0	08 - 1,46	gw / borderline		
> 20,0	> 0,88	> 0,92	> 1,00	> 1,08	> 1,15	> 1,23	> 1,31	> 1,38	> 1,46	pos		

Formeln für spezielle Auswertesysteme Special case formulas OD = 1,023 x MV(STD) entspricht oberem cut-off/corresponds to upper cut-off OD = 0,754 x MV(STD) entspricht unterem cut-off/corresponds to lower cut-off

 $Concentration = \exp(3,881 - \ln(3,693/(\text{MV(Sample}) \ x0,86/\ \text{MV(STD)} - 0,09) - 1)/1,471)$

20 Institut Virion\Serion GmbH15 Friedrich-Bergius-Ring 19

Friedrich-Bergius-Ring 19 D-97076 Würzburg

Zusätzliche Barcodes mit Formeln für / Additional Barcodes with formulas for Revelation™ DSX / DS-Matrix ™

4PS- Formel / 4PS-formula

exp(3.881-ln(3.693/(Sample*0.860/S-0.090)-1)/1.471)



Gültigkeitsbereich / Validity Range

0.430<=S1<=1.462



If OD Sample < Parameter A

if Ti<(0.090*(S1/0.860)) then Ti=(0.090+0.001)*(S1/0.860)



If OD Sample > Parameter D

if Ti>(3.783*(S1/0.860)) then Ti=(3.783-0.001)*(S1/0.860)



If OD Negative control < Parameter A

if NC1 < (0.090*(S1/0.860)) then NCi = (0.090+0.001)*(S1/0.860)



Institut Virion\Serion GmbH Friedrich-Bergius-Ring 19 D-97076 Würzburg