SERION ELISA classic ESR112G

FSME VIRUS / TBE VIRUS IgG

SMH.BI

Qualitätskontrollzertifikat / Quality Control Certificate

Kitcharge / Lot SMH.BI IFU-Version 112-15

Verw. bis / Exp. 2019-11

22.12.2017

Prüfdatum /

Date of control



150

100

Verwendete Reagenzien / Reagents used	Standard							Standard Kur	Standard Kurve / Standard curve		
Teststreifen / Antigen coated strips	SIH.CI	Ref Werte / Re	f. Values		Gültigke	itsbereich	/ Valid	lity Range	Parameter	Α	-0,011
Standardserum / Standard serum	SLH.DW	10	0,86		OD	0,43	-	1,46		В	0,980
Negativ Kontrolle / Negative control	SLH.DV									С	6,498
Konjugat / Conjugate	SHH.EA+++	Units	s 524 U /m	l						D	1,959
Quantifizierungsgrenzen / Limits of quantific	U/ml	30		3000							
Grenzwertbereich / Borderline range	U/ml	100	-	150							

OD Bereich / OD Range 405 nm, Standardserum / Standard serum																			
0,43	} -	0,47	0,48	- 0,53	0,54	- 0,58	0,59	- 0,64	0,65	- 0,69	0,70	- 0,74	0,75	- 0,80	0,81	- 0,85	0,86	U/ml	Interpretation
	<	0,14		< 0,15		< 0,17		< 0,19		< 0,20		< 0,22		< 0,23		< 0,25	< 0,26	< 100,0	neg
0,14		0,19	0,15	- 0,21	0,17	- 0,23	0,19	- 0,26	0,20	- 0,28	0,22	- 0,30	0,23	- 0,32	0,25	- 0,35	0,26 - 0,36	100,0 - 150,0	gw / borderline
	>	0,19		> 0,21		> 0,23		> 0,26		> 0,28		> 0,30		> 0,32		> 0,35	> 0,36	> 150,0	pos

OD Bereich / OD Range 405 nm, Standardserum / Standard serum											
U/ml	0,86	0,87 - 0,94	0,95 - 1,01	1,02 - 1,09	1,10 - 1,16	1,17 - 1,24	1,25 - 1,31	1,32 - 1,39 1,40 - 1,46	Interpretation		
< 100,0	< 0,26	< 0,27	< 0,30	< 0,32	< 0,34	< 0,36	< 0,39	< 0,41 < 0,43	neg		
100,0 - 150,0 > 150,0	, ,	0,27 - 0,38 > 0,38	0,30 - 0,41 > 0,41	0,32 - 0,44 > 0,44	0,34 - 0,47 > 0,47	0,36 - 0,50 > 0,50	0,39 - 0,54 > 0,54		gw / borderline pos		

Formeln für spezielle Auswertesysteme Special case formulas OD = **0,420** x MV(STD) entspricht oberem cut-off/corresponds to upper cut-off
OD = **0,297** x MV(STD) entspricht unterem cut-off/corresponds to lower cut-off

 $Concentration = exp(6,498\text{-In}(1,97/(MV(Sample)\ x0,86/\ MV(STD) + 0,011) - 1)/0,98)$

Institut Virion\Serion GmbH Friedrich-Bergius-Ring 19 D-97076 Würzburg

Zusätzliche Barcodes mit Formeln für / Additional Barcodes with formulas for Revelation™ DSX / DS-Matrix ™

4PS- Formel / 4PS-formula

exp(6.498-ln(1.970/(Sample*0.860/S+0.011)-1)/0.980)



Gültigkeitsbereich / Validity Range

0.430<=S1<=1.462



If OD Sample < Parameter A

if Ti < (-0.011*(S1/0.860)) then Ti = (-0.011+0.001)*(S1/0.860)



If OD Sample > Parameter D

if Ti>(1.959*(S1/0.860)) then Ti=(1.959-0.001)*(S1/0.860)



If OD Negative control < Parameter A

if NC1<(-0.011*(S1/0.860)) then NCi=(-0.011+0.001)*(S1/0.860)



Institut Virion\Serion GmbH Friedrich-Bergius-Ring 19 D-97076 Würzburg