SERION ELISA classic

ESR1361M

EPSTEIN-BARR VIRUS / VCA IgM

SGI.DK

Qualitätskontrollzertifikat / Quality Control Certificate

Kitcharge / Lot SGI.DK IFU-Version 136-23

Verw. bis / Exp. 2020-03

31.07.2018

Prüfdatum /

Date of control



13

9

Verwendete Reagenzien / Reagents used	Standard		Standard Kur	Standard Kurve / Standard curve							
Teststreifen / Antigen coated strips	Ref Werte / Ref. Values			Gültigkeitsbereich / Validity Range				Parameter	Α	-0,032	
Standardserum / Standard serum	SDI.AW	OD	0,96		OD	0,48	-	1,63		В	0,770
Negativ Kontrolle / Negative control	SDI.AV									С	5,307
Konjugat / Conjugate	SEI.BS++	Units	34,4	U/ml						D	4,835
Quantifizierungsgrenzen / Limits of quantific	U/ml	4	-	200							
Grenzwertbereich / Borderline range	U/ml	9	-	13							

OD Bereich / OD Range 405 nm, Standardserum / Standard serum																				
0,48	-	0,53	0,54	- 0,59	0,60	- 0,65	0,66	- 0,71	0,72	- 0,77	0,78	- 0,83	0,84	- 0,89	0,90	- 0,95	0,96	U/ml		Interpretation
	<	0,19		< 0,22		< 0,24		< 0,26		< 0,29		< 0,31		< 0,33		< 0,36	< 0,37	<	9,0	neg
0,19	-	0,26	0,22	- 0,29	0,24	- 0,32	0,26	- 0,35	0,29	- 0,38	0,31	- 0,41	0,33	- 0,44	0,36	- 0,47	0,37 - 0,49	9,0 -	13,0	gw / borderline
	>	0,26		> 0,29		> 0,32		> 0,35		> 0,38		> 0,41		> 0,44		> 0,47	> 0,49	>	13,0	pos

OD Bereich / OD Range 405 nm, Standardserum / Standard serum											
U/ml			0,96	0,97 - 1,04	1,05 - 1,13	1,14 - 1,21	1,22 - 1,30	1,31 - 1,38	1,39 - 1,46	1,47 - 1,55 1,56 - 1,63	Interpretation
	<	9,0	< 0,37	< 0,39	< 0,42	< 0,45	< 0,49	< 0,52	< 0,55	< 0,58 < 0,61	neg
9,0	-	13,0 13,0	, ,	0,39 - 0,51 > 0,51	0,42 - 0,56 > 0,56	0,45 - 0,60 > 0,60	0,49 - 0,64 > 0,64	0,52 - 0,69 > 0,69	0,55 - 0,73 > 0,73		gw / borderline pos

Formeln für spezielle Auswertesysteme Special case formulas $\begin{array}{lll} \text{OD} = & \textbf{0,514} & \text{x MV(STD)} \text{ entspricht oberem cut-off/corresponds to upper cut-off} \\ \text{OD} = & \textbf{0,390} & \text{x MV(STD)} \text{ entspricht unterem cut-off/corresponds to lower cut-off} \\ \end{array}$

 $Concentration = \exp(5,307 - \ln(4,867/(\text{MV(Sample}) \ x0,96/\ \text{MV(STD)} + 0,032) - 1)/0,77)$

Institut Virion\Serion GmbH

Friedrich-Bergius-Ring 19 D-97076 Würzburg

Zusätzliche Barcodes mit Formeln für / Additional Barcodes with formulas for Revelation™ DSX / DS-Matrix ™

4PS- Formel / 4PS-formula

exp(5.307-ln(4.867/(Sample*0.960/S+0.032)-1)/0.770)



Gültigkeitsbereich / Validity Range

0.480<=S1<=1.632



If OD Sample < Parameter A

if Ti < (-0.032*(S1/0.960)) then Ti = (-0.032+0.001)*(S1/0.960)



If OD Sample > Parameter D

if Ti>(4.835*(S1/0.960)) then Ti=(4.835-0.001)*(S1/0.960)



If OD Negative control < Parameter A

if NC1<(-0.032*(S1/0.960)) then NCi=(-0.032+0.001)*(S1/0.960)



Institut Virion\Serion GmbH Friedrich-Bergius-Ring 19 D-97076 Würzburg