SERION ELISA classic ESR134A

COXSACKIEVIRUS IgA

SCI.BU

Qualitätskontrollzertifikat / Quality Control Certificate

Kitcharge / Lot SCI.BU IFU-Version 134-15

Verw. bis / Exp. 2020-02

16.03.2018

Prüfdatum /

Date of control



15

10

Verwendete Reagenzien / Reagents used	Standard								Standard Ku	Standard Kurve / Standard curve		
Teststreifen / Antigen coated strips	SBH.CV	Ref Werte / Re	f. Value	es	Gü	Itigkei	tsbereich	/ Valid	dity Range	Parameter	Α	0,009
Standardserum / Standard serum	SBI.EI	10	0,92			DD	0,46	-	1,56		В	1,082
Negativ Kontrolle / Negative control	SBI.EH										С	4,541
Konjugat / Conjugate	SIH.CE+	Unit	s 24,1	U/ml							D	4,886
Quantifizierungsgrenzen / Limits of quantific	U/ml	4	-		150							
Grenzwertbereich / Borderline range	U/ml	10) -		15							

OD Bereich / OD Range 405 nm, Standardserum / Standard serum																			
0,46	-	0,51	0,52	- 0,57	0,58	- 0,62	0,63	- 0,68	0,69	- 0,74	0,75	- 0,80	0,81	- 0,85	0,86	- 0,91	0,92	U/ml	Interpretation
	<	0,22		< 0,24		< 0.27		< 0.29		< 0.32		< 0.34		< 0.37		< 0.39	< 0.41	< 10	.0 neg
0,22	-		0,24		0,27	- 0,39	0,29	,	0,32	- 0,47	0,34	,	0,37	- 0,54	0,39	- 0,58	0,41 - 0,60		,0 gw / borderline
	>	0,32		> 0,35		> 0,39		> 0,43		> 0,47		> 0,50		> 0,54		> 0,58	> 0,60		,0 pos

OD Bereich / OD Range 405 nm, Standardserum / Standard serum										
U/ml 0,92	0,93 - 1,00	1,01 - 1,08	1,09 - 1,16	1,17 - 1,24	1,25 - 1,32	1,33 - 1,40	1,41 - 1,48 1,49 - 1,56	Interpretation		
< 10,0 <	0,41 < 0,43	< 0,47	< 0,50	< 0,54	< 0,57	< 0,61	< 0,65 < 0,68	neg		
10,0 - 15,0 0,41 - > 15.0 >	0,60	0,47 - 0,68 > 0,68	0,50 - 0,73 > 0,73	0,54 - 0,79 > 0,79	0,57 - 0,84 > 0,84	0,61 - 0,89 > 0.89		gw / borderline		

Formeln für spezielle Auswertesysteme Special case formulas OD = **0,651** x MV(STD) entspricht oberem cut-off/corresponds to upper cut-off
OD = **0,442** x MV(STD) entspricht unterem cut-off/corresponds to lower cut-off

Concentration= exp(4,541-In(4,877/(MV(Sample) x0,92/ MV(STD)-0,009)-1)/1,082)

Institut Virion\Serion GmbH

Friedrich-Bergius-Ring 19 D-97076 Würzburg

Zusätzliche Barcodes mit Formeln für / Additional Barcodes with formulas for Revelation™ DSX / DS-Matrix ™

4PS- Formel / 4PS-formula

exp(4.541-ln(4.877/(Sample*0.920/S-0.009)-1)/1.082)



Gültigkeitsbereich / Validity Range

0.460<=S1<=1.564



If OD Sample < Parameter A

if Ti < (0.009*(S1/0.920)) then Ti = (0.009+0.001)*(S1/0.920)



If OD Sample > Parameter D

if Ti>(4.886*(S1/0.920)) then Ti=(4.886-0.001)*(S1/0.920)



If OD Negative control < Parameter A

if NC1<(0.009*(S1/0.920)) then NCi=(0.009+0.001)*(S1/0.920)



Institut Virion\Serion GmbH Friedrich-Bergius-Ring 19 D-97076 Würzburg