## **SERION ELISA classic**

**ESR1231M** 

## **INFLUENZA A VIRUS IgM**

SBI.DS

## **Qualitätskontrollzertifikat / Quality Control Certificate**

Kitcharge / Lot SBI.DS IFU-Version 123-15

20.02.2018 Prüfdatum /

Verw. bis / Exp. 2020-02

Date of control



15

10

Verwendete Reagenzien / Reagents used	Standard		Standard Kur	Standard Kurve / Standard curve						
Teststreifen / Antigen coated strips	Ref Werte / Re	f. Values		Gültigkei	tsbereich	/ Validity Range	Parameter	Α	0,067	
Standardserum / Standard serum	OD <b>0,86</b>			OD <b>0,43</b> - <b>1,46</b>				В	1,220	
Negativ Kontrolle / Negative control	SBI.AN								С	4,360
Konjugat / Conjugate	SAI.DN+++	Units	s <b>23,2 U</b> /	ml					D	4,351
Quantifizierungsgrenzen / Limits of quantific	U/ml	5	-	100						
Grenzwertbereich / Borderline range	U/ml	10	-	15						

OD Bereich / OD Range 405 nm, Standardserum / Standard serum																					
0,43	-	0,47	0,48	- 0,53	0,54	- 0,58	0,59	- 0,64	0,65	- 0,69	0,70	- 0,74	0,75	- 0,80	0,81	- 0,85	0,86	U/	/ml		Interpretation
	<	0,20		< 0.23		< 0.25		< 0.28		< 0.30		< 0.33		< 0.35		< 0.38	< 0.	39	<	10,0	neg
0,20	-	0,30	0,23	- 0,34	0,25	- 0,37	0,28	- 0,41	0,30	- 0,44	0,33	- 0,48	0,35	- 0,51	0,38	- 0,55	0,39 - 0,	57	10,0 -	15,0	gw / borderline
	>	0,30		> 0,34		> 0,37		> 0,41		> 0,44		> 0,48		> 0,51		> 0,55	> 0.	57	>	15,0	pos

OD Bereich / OD Range 405 nm, Standardserum / Standard serum												
U/ml		0,86	0,87 - 0,94	0,95 - 1,01	1,02 - 1,09	1,10 - 1,16	1,17 - 1,24	1,25 - 1,31	1,32 - 1,39 1,40 - 1,46	Interpretation		
<	10,0	< 0,39	< 0,41	< 0,44	< 0,48	< 0,51	< 0,55	< 0,58	< 0,61 < 0,65	neg		
10,0 -	15,0 15,0	, ,	0,41 - 0,60 > 0,60	0,44 - 0,65 > 0,65	0,48 - 0,70 > 0,70	0,51 - 0,75 > 0,75	0,55 - 0,80 > 0,80	0,58 - 0,85 > 0,85		gw / borderline pos		

Formeln für spezielle Auswertesysteme Special case formulas OD = **0,664** x MV(STD) entspricht oberem cut-off/corresponds to upper cut-off
OD = **0,452** x MV(STD) entspricht unterem cut-off/corresponds to lower cut-off

 $Concentration = \exp(4,36 - \ln(4,284 / (MV(Sample) \ x0,86 / \ MV(STD) - 0,067) - 1) / 1,22)$ 

Institut Virion\Serion GmbH

Friedrich-Bergius-Ring 19 D-97076 Würzburg

## Zusätzliche Barcodes mit Formeln für / Additional Barcodes with formulas for Revelation™ DSX / DS-Matrix ™

4PS- Formel / 4PS-formula

exp(4.360-ln(4.284/(Sample\*0.860/S-0.067)-1)/1.220)



Gültigkeitsbereich / Validity Range

0.430<=S1<=1.462



If OD Sample < Parameter A

if Ti < (0.067\*(S1/0.860)) then Ti = (0.067+0.001)\*(S1/0.860)



If OD Sample > Parameter D

if Ti>(4.351\*(S1/0.860)) then Ti=(4.351-0.001)\*(S1/0.860)



If OD Negative control < Parameter A

if NC1<(0.067\*(S1/0.860)) then NCi=(0.067+0.001)\*(S1/0.860)



Institut Virion\Serion GmbH Friedrich-Bergius-Ring 19 D-97076 Würzburg