SERION ELISA classic

ESR129M

ROETELN VIRUS / RUBELLA VIRUS IgM

14.12.2017

SMH.BB

Qualitätskontrollzertifikat / Quality Control Certificate

Kitcharge / Lot SMH.BB IFU-Version 129-16

Verw. bis / Exp.

Prüfdatum / 2019-11

Date of control



3,5

2,5

Verwendete Reagenzien / Reagents used	Standard		Standard Kur	Standard Kurve / Standard curve							
Teststreifen / Antigen coated strips	SLH.AR	Ref Werte / Re	f. Values		Gültigke	itsbereich	/ Valid	dity Range	Parameter	Α	0,018
Standardserum / Standard serum	SLH.DU	OD	0,81		OD	0,41	-	1,38		В	0,922
Negativ Kontrolle / Negative control	SLH.DT									С	3,843
Konjugat / Conjugate	SKH.AQ++	Units	9,65 U/	ml						D	4,200
Quantifizierungsgrenzen / Limits of quantific	U/ml	2	-	120							
Grenzwertbereich / Borderline range	U/ml	2,5	-	3,5							

OD Bereich / OD Range 405 nm, Standardserum / Standard serum																								
0,41		- (0,45	0,46	- 0,50	0,51	- 0,55	0,56	- 0,60	0,61	- 0	0,65	0,66	- 0,70	0,71	- 0,75	0,76	- 0	,80	0,81	U/ml			Interpretation
	<	< (0,15		< 0,16		< 0,18		< 0,20		< 0	0,22		< 0,23		< 0,25		< 0	,27	< 0,28		<	2,5	neg
0,15	5 .	- (0,19	0,16	- 0,22	0,18	- 0,24	0,20	- 0,26	0,22	- 0	0,29	0,23	- 0,31	0,25	- 0,33	0,27	- 0	,36	0,28 - 0,37	2,5	-	3,5	gw / borderline
	>	> (0,19		> 0,22		> 0,24		> 0,26		> 0	0,29		> 0,31		> 0,33		> 0	,36	> 0,37		>	3,5	pos

1,25 - 1,31 1,32 - 1,38	Interpretation
< 0,44 < 0,47	neg
	gw / borderline
	< 0,44 < 0,47

Formeln für spezielle Auswertesysteme Special case formulas

0,456 x MV(STD) entspricht oberem cut-off/corresponds to upper cut-off OD = 0,348 x MV(STD) entspricht unterem cut-off/corresponds to lower cut-off OD =

Concentration= exp(3,843-In(4,182/(MV(Sample) x0,81/ MV(STD)-0,018)-1)/0,922)

Institut Virion\Serion GmbH Friedrich-Bergius-Ring 19 D-97076 Würzburg

Zusätzliche Barcodes mit Formeln für / Additional Barcodes with formulas for Revelation™ DSX / DS-Matrix ™

4PS- Formel / 4PS-formula

exp(3.843-ln(4.182/(Sample*0.810/S-0.018)-1)/0.922)



Gültigkeitsbereich / Validity Range

0.405<=S1<=1.377



If OD Sample < Parameter A

if Ti < (0.018*(S1/0.810)) then Ti = (0.018+0.001)*(S1/0.810)



If OD Sample > Parameter D

if Ti>(4.200*(S1/0.810)) then Ti=(4.200-0.001)*(S1/0.810)



If OD Negative control < Parameter A

if NC1 < (0.018*(S1/0.810)) then NCi = (0.018+0.001)*(S1/0.810)



Institut Virion\Serion GmbH Friedrich-Bergius-Ring 19 D-97076 Würzburg