SERION ELISA classic

ESR1231G

INFLUENZA A VIRUS IgG Kinder/Children

SBI.EZ

Qualitätskontrollzertifikat / Quality Control Certificate

Kitcharge / Lot SBI.EZ IFU-Version 123-15

2020-02

Verw. bis / Exp.

02.03.2018 Prüfdatum /

Date of control



Verwendete Reagenzien / Reagents used	Standard			Standard Kur	Standard Kurve / Standard curve						
Teststreifen / Antigen coated strips	SCH.AI	Ref Werte / Ref	f. Values		Gültigke	itsbereich	/ Valid	dity Range	Parameter	Α	0,010
Standardserum / Standard serum	SBI.CI	OD	0,78		OD	0,39	-	1,33		В	1,018
Negativ Kontrolle / Negative control	SBI.CH									С	3,420
Konjugat / Conjugate	SLH.CD++	Units	30,9 l	J/ml						D	1,543
Quantifizierungsgrenzen / Limits of quantific	U/ml	2	-	200							
Grenzwertbereich / Borderline range	U/ml	4	-	7							

OD Bereich / OD Range 405 nm, Standardserum / Standard serum																			
0,39	-	0,43	0,44	- 0,48	0,49	- 0,53	0,54	- 0,58	0,59	- 0,62	0,63	- 0,67	0,68	- 0,72	0,73	- 0,77	0,78	U/ml	Interpretation
	<	0,09		< 0,11		< 0,12		< 0,13		< 0,14		< 0,15		< 0,16		< 0,17	< 0,18	< 4,0	neg
0,09	-	0,15	0,11	- 0,17	0,12	- 0,19	0,13	- 0,21	0,14	- 0,22	0,15	- 0,24	0,16	- 0,26	0,17	- 0,28	0,18 - 0,29	4,0 - 7,0	gw / borderline
	>	0,15		> 0,17		> 0,19		> 0,21		> 0,22		> 0,24		> 0,26		> 0,28	> 0,29	> 7,0	pos

11/11/1 0.70 0.05 0.00 0.00 0.00 0.00 4.05 4.00 4.40 4.4		
U/ml 0,78 0,79 - 0,85 0,86 - 0,92 0,93 - 0,98 0,99 - 1,05 1,06 - 1,12 1,13 - 1,19	1,20 - 1,26 1,27 - 1,33	Interpretation
< 4,0 < 0,18 < 0,19 < 0,20 < 0,22 < 0,24 < 0,25 < 0,27	< 0,28 < 0,30	neg
4,0 - 7,0 0,18 - 0,29 0,19 - 0,30 0,20 - 0,33 0,22 - 0,36 0,24 - 0,38 0,25 - 0,41 0,27 - 0,43 > 7,0 > 0,29 > 0,30 > 0,33 > 0,36 > 0,38 > 0,41 > 0,43		gw / borderline pos

Formeln für spezielle Auswertesysteme Special case formulas

0,371 x MV(STD) entspricht oberem cut-off/corresponds to upper cut-off OD =

0,233 x MV(STD) entspricht unterem cut-off/corresponds to lower cut-off OD =

Concentration= exp(3,42-ln(1,533/(MV(Sample) x0,78/ MV(STD)-0,01)-1)/1,018)

Institut Virion\Serion GmbH 7

4

Friedrich-Bergius-Ring 19 D-97076 Würzburg

Zusätzliche Barcodes mit Formeln für / Additional Barcodes with formulas for Revelation™ DSX / DS-Matrix ™

4PS- Formel / 4PS-formula

exp(3.420-ln(1.533/(Sample*0.780/S-0.010)-1)/1.018)



Gültigkeitsbereich / Validity Range

0.390<=S1<=1.326



If OD Sample < Parameter A

if Ti < (0.010*(S1/0.780)) then Ti = (0.010+0.001)*(S1/0.780)



If OD Sample > Parameter D

if Ti>(1.543*(S1/0.780)) then Ti=(1.543-0.001)*(S1/0.780)



If OD Negative control < Parameter A

if NC1 < (0.010*(S1/0.780)) then NCi = (0.010+0.001)*(S1/0.780)



Institut Virion\Serion GmbH Friedrich-Bergius-Ring 19 D-97076 Würzburg