SERION ELISA classic

ESR1361G

EPSTEIN-BARR VIRUS / VCA IgG

SGI.ED

Qualitätskontrollzertifikat / Quality Control Certificate

Kitcharge / Lot SGI.ED IFU-Version 136-23

Verw. bis / Exp. 2020-07

30.07.2018

Prüfdatum /

Date of control



Verwendete Reagenzien / Reagents used	Lot	Standard	tandard							Standard Kurve / Standard curve		
Teststreifen / Antigen coated strips	SGI.AS	Ref Werte / Re	f. Values		Gültigke	itsbereich	/ Valid	lity Range	Parameter	Α	-0,019	
Standardserum / Standard serum	SGI.CM	OE	0,85		OD	0,43	-	1,45		В	0,842	
Negativ Kontrolle / Negative control	SGI.CI									С	4,745	
Konjugat / Conjugate	SDI.EB+++	Units	s 64,5 U	J/ml						D	2,265	
Quantifizierungsgrenzen / Limits of quantific	U/ml	4	-	200								
Grenzwertbereich / Borderline range	U/ml	10	-	15								

OD Bereich / OD Range 405 nm, Standardserum / Standard serum																			
0,43	-	0,47	0,48	- 0,52	0,53	- 0,57	0,58	- 0,63	0,64	- 0,68	0,69	- 0,73	0,74	- 0,79	0,80	- 0,84	0,85	U/ml	Interpretation
	<	0,13		< 0.14		< 0.16		< 0,17		< 0.19		< 0.20		< 0,22		< 0.23	< 0,24	< 10.0	neg
0,13	-	0,17	0,14	- 0,19	0,16	- 0,21	0,17	- 0,24	0,19	- 0,26	0,20	- 0,28	0,22	- 0,30	0,23	- 0,32	0,24 - 0,33	10,0 - 15,0	gw / borderline
	>	0,17		> 0,19		> 0,21		> 0,24		> 0,26		> 0,28		> 0,30		> 0,32	> 0,33	> 15,0	pos

OD Bereich / OD Range 405 nm, Standardserum / Standard serum												
U/ml	0,85	0,86 - 0,92	0,93 - 1,00	1,01 - 1,07	1,08 - 1,15	1,16 - 1,22	1,23 - 1,30	1,31 - 1,37 1,38 - 1,45	Interpretation			
< 10,0	< 0,24	< 0,25	< 0,27	< 0,29	< 0,31	< 0,34	< 0,36	< 0,38 < 0,40	neg			
10,0 - 15,0 > 15,0	0,24 - 0,33 > 0,33	0,25 - 0,35 > 0,35	0,27 - 0,38 > 0,38	0,29 - 0,40 > 0,40	0,31 - 0,43 > 0,43	0,34 - 0,46 > 0,46	0,36 - 0,49 > 0.49		gw / borderline pos			

Formeln für spezielle Auswertesysteme Special case formulas OD = **0,387** x MV(STD) entspricht oberem cut-off/corresponds to upper cut-off
OD = **0,282** x MV(STD) entspricht unterem cut-off/corresponds to lower cut-off

Concentration= exp(4,745-ln(2,284/(MV(Sample) x0,85/ MV(STD)+0,019)-1)/0,842)

15 Institut Virion\Serion GmbH

10

Friedrich-Bergius-Ring 19 D-97076 Würzburg

Zusätzliche Barcodes mit Formeln für / Additional Barcodes with formulas for Revelation™ DSX / DS-Matrix ™

4PS- Formel / 4PS-formula

exp(4.745-ln(2.284/(Sample*0.850/S+0.019)-1)/0.842)



Gültigkeitsbereich / Validity Range

0.425<=S1<=1.445



If OD Sample < Parameter A

if Ti < (-0.019*(S1/0.850)) then Ti = (-0.019+0.001)*(S1/0.850)



If OD Sample > Parameter D

if Ti>(2.265*(S1/0.850)) then Ti=(2.265-0.001)*(S1/0.850)



If OD Negative control < Parameter A

if NC1<(-0.019*(S1/0.850)) then NCi=(-0.019+0.001)*(S1/0.850)



Institut Virion\Serion GmbH Friedrich-Bergius-Ring 19 D-97076 Würzburg