SERION ELISA classic ESR118M

HELICOBACTER PYLORI IgM

SKH.DD

Qualitätskontrollzertifikat / Quality Control Certificate

Kitcharge / Lot SKH.DD IFU-Version 118-14

Verw. bis / Exp. 09-2019 Prüfdatum /

Date of control

02.11.2017



30

20

Verwendete Reagenzien / Reagents used	Standard							Standard Kur	Standard Kurve / Standard curve		
Teststreifen / Antigen coated strips	Ref Werte / Ref	f. Values		Gültigke	itsbereich	/ Valid	dity Range	Parameter	Α	0,003	
Standardserum / Standard serum	OD	0,93		OD	0,47	-	1,58		В	0,969	
Negativ Kontrolle / Negative control	SKH.CB									С	5,263
Konjugat / Conjugate	SHH.BD+	Units	68,8 U/ı	nl						D	3,451
Quantifizierungsgrenzen / Limits of quantific	U/ml	10	-	300							
Grenzwertbereich / Borderline range	U/ml	20	-	30							

OD Bereich / OD Range 405 nm, Standardserum / Standard serum																			
0,47		0,51	0,52	- 0,57	0,58	- 0,63	0,64	- 0,69	0,70	- 0,75	0,76	- 0,80	0,81	- 0,86	0,87	- 0,92	0,93	U/ml	Interpretation
	<	: 0,18		< 0.21		< 0.23		< 0.25		< 0.27		< 0.29		< 0.32		< 0.34	< 0.35	< 20.0	neg
0,18	} -	0,26	0,21	- 0,29	0,23	- 0,32	0,25	- 0,35	0,27	- 0,38	0,29	,	0,32	- 0,44	0,34	- 0,47	0,35 - 0,49		gw / borderline
	>	> 0,26		> 0,29		> 0,32		> 0,35		> 0,38		> 0,41		> 0,44		> 0,47	> 0,49	> 30,0	pos

OD Bereich / OD Range 405 nm, Standardserum / Standard serum												
U/ml		0,93	0,94 - 1,01	1,02 - 1,09	1,10 - 1,17	1,18 - 1,26	1,27 - 1,34	1,35 - 1,42	1,43 - 1,50 1,51 - 1,58	Interpretation		
	< 20,0	< 0,35	< 0,37	< 0,40	< 0,43	< 0,46	< 0,49	< 0,52	< 0,55 < 0,58	neg		
20,0	- 30,0 > 30,0	, ,	0,37 - 0,51 > 0,51	0,40 - 0,56 > 0,56	0,43 - 0,60 > 0,60	0,46 - 0,64 > 0,64	0,49 - 0,69 > 0,69	0,52 - 0,73 > 0,73	0,55 - 0,77 0,58 - 0,81 > 0,77 > 0,81	gw / borderline pos		

Formeln für spezielle Auswertesysteme Special case formulas OD = **0,527** x MV(STD) entspricht oberem cut-off/corresponds to upper cut-off
OD = **0,374** x MV(STD) entspricht unterem cut-off/corresponds to lower cut-off

 $Concentration = exp(5,263-ln(3,448/(MV(Sample)\ x0,93/\ MV(STD)-0,003)-1)/0,969)$

Institut Virion\Serion GmbH Friedrich-Bergius-Ring 19 D-97076 Würzburg

Zusätzliche Barcodes mit Formeln für / Additional Barcodes with formulas for Revelation™ DSX / DS-Matrix ™

4PS- Formel / 4PS-formula

exp(5.263-ln(3.448/(Sample*0.930/S-0.003)-1)/0.969)



Gültigkeitsbereich / Validity Range

0.465<=S1<=1.581



If OD Sample < Parameter A

if Ti<(0.003*(S1/0.930)) then Ti=(0.003+0.001)*(S1/0.930)



If OD Sample > Parameter D

if Ti>(3.451*(S1/0.930)) then Ti=(3.451-0.001)*(S1/0.930)



If OD Negative control < Parameter A

if NC1<(0.003*(S1/0.930)) then NCi=(0.003+0.001)*(S1/0.930)



Institut Virion\Serion GmbH Friedrich-Bergius-Ring 19 D-97076 Würzburg