**SERION ELISA classic** 

**ESR1363G** 

## **EPSTEIN-BARR VIRUS / EA IgG**

SCI.BB

### **Qualitätskontrollzertifikat / Quality Control Certificate**

Kitcharge / Lot

Verw. bis / Exp.

SCI.BB

2019-11

IFU-Version 136-23

06.03.2018

Prüfdatum /

Date of control



15

10

| Verwendete Reagenzien / Reagents used         | Standard |                 |          |      |          |             |        | Standard Ku | Standard Kurve / Standard curve |   |        |
|---|----------|-----------------|----------|------|----------|-------------|--------|-------------|---------------------------------|---|--------|
| Teststreifen / Antigen coated strips          | SLH.DG   | Ref Werte / Ref | f. Value | s    | Gültigke | eitsbereich | / Vali | idity Range | Parameter                       | Α | -0,019 |
| Standardserum / Standard serum                | SAI.FC   | OD              | 0,81     |      | OD       | 0,41        | -      | 1,38        |                                 | В | 0,958  |
| Negativ Kontrolle / Negative control          | SAI.FB   |                 |          |      |          |             |        |             |                                 | С | 4,405  |
| Konjugat / Conjugate                          | SLH.CD++ | Units           | 69,6     | U/ml |          |             |        |             |                                 | D | 1,778  |
| Quantifizierungsgrenzen / Limits of quantific | U/ml     | 2               | -        | 400  |          |             |        |             |                                 |   |        |
| Grenzwertbereich / Borderline range           | U/ml     | 10              | -        | 15   |          |             |        |             |                                 |   |        |

| OD Bereich / OD Range 405 nm, Standardserum / Standard serum |     |      |      |        |      |        |      |        |      |        |      |        |      |        |      |        |             |             |                 |
|--|-----|------|------|--------|------|--------|------|--------|------|--------|------|--------|------|--------|------|--------|-------------|-------------|-----------------|
| 0,41   | -   | 0,45 | 0,46 | - 0,50 | 0,51 | - 0,55 | 0,56 | - 0,60 | 0,61 | - 0,65 | 0,66 | - 0,70 | 0,71 | - 0,75 | 0,76 | - 0,80 | 0,81        | U/ml        | Interpretation  |
|  | <   | 0,10 |      | < 0,11 |      | < 0,12 |      | < 0,14 |      | < 0,15 |      | < 0,16 |      | < 0,17 |      | < 0,18 | < 0,19      | < 10,0      | neg             |
| 0,10   | ) - | 0,15 | 0,11 | - 0,16 | 0,12 | - 0,18 | 0,14 | - 0,20 | 0,15 | - 0,22 | 0,16 | - 0,23 | 0,17 | - 0,25 | 0,18 | - 0,27 | 0,19 - 0,28 | 10,0 - 15,0 | gw / borderline |
|  | >   | 0,15 |      | > 0,16 |      | > 0,18 |      | > 0,20 |      | > 0,22 |      | > 0,23 |      | > 0,25 |      | > 0,27 | > 0,28      | > 15,0      | pos             |

| OD Bereich / OD Range 405 nm, Standardserum / Standard serum |                       |                       |                       |                       |                       |                       |                       |                         |                 |  |  |
|--|-----------------------|-----------------------|-----------------------|-----------------------|-----------------------|-----------------------|-----------------------|-------------------------|-----------------|--|--|
| U/ml   | 0,81                  | 0,82 - 0,88           | 0,89 - 0,95           | 0,96 - 1,02           | 1,03 - 1,09           | 1,10 - 1,16           | 1,17 - 1,24           | 1,25 - 1,31 1,32 - 1,38 | Interpretation  |  |  |
| < 10,0   | < 0,19                | < 0,20                | < 0,22                | < 0,23                | < 0,25                | < 0,27                | < 0,28                | < 0,30 < 0,32           | neg             |  |  |
| 10,0 - 15,0<br>> 15,0  | 0,19 - 0,28<br>> 0,28 | 0,20 - 0,29<br>> 0,29 | 0,22 - 0,32<br>> 0,32 | 0,23 - 0,34<br>> 0,34 | 0,25 - 0,37<br>> 0,37 | 0,27 - 0,39<br>> 0,39 | 0,28 - 0,42<br>> 0.42 |                         | gw / borderline |  |  |

Formeln für spezielle Auswertesysteme Special case formulas  $\mathsf{OD} = \qquad \textbf{0,341} \qquad \mathsf{x} \; \mathsf{MV}(\mathsf{STD}) \; \mathsf{entspricht} \; \mathsf{oberem} \; \mathsf{cut\text{-}off/corresponds} \; \mathsf{to} \; \mathsf{upper} \; \mathsf{cut\text{-}off}$ 

OD =  $0.238 \times MV(STD)$  entspricht unterem cut-off/corresponds to lower cut-off Concentration= exp(4.405-ln(1.797/(MV(Sample) x0.81/MV(STD)+0.019)-1)/0.958)

Institut Virion\Serion GmbH

Friedrich-Bergius-Ring 19 D-97076 Würzburg

# Zusätzliche Barcodes mit Formeln für / Additional Barcodes with formulas for Revelation™ DSX / DS-Matrix ™

#### 4PS- Formel / 4PS-formula

exp(4.405-ln(1.797/(Sample\*0.810/S+0.019)-1)/0.958)



#### Gültigkeitsbereich / Validity Range

0.405<=S1<=1.377



#### If OD Sample < Parameter A

if Ti < (-0.019\*(S1/0.810)) then Ti = (-0.019+0.001)\*(S1/0.810)



#### If OD Sample > Parameter D

if Ti>(1.778\*(S1/0.810)) then Ti=(1.778-0.001)\*(S1/0.810)



#### If OD Negative control < Parameter A

if NC1 < (-0.019\*(S1/0.810)) then NCi = (-0.019+0.001)\*(S1/0.810)



Institut Virion\Serion GmbH Friedrich-Bergius-Ring 19 D-97076 Würzburg