SERION ELISA classic

ESR1231A

INFLUENZA A VIRUS IgA Erwachsene/Adults

SCI.BA

Qualitätskontrollzertifikat / Quality Control Certificate

Kitcharge / Lot SCI.BA IFU-Version 123-15

Verw. bis / Exp. 2020-03

05.03.2018

Prüfdatum /

Date of control



15

10

Verwendete Reagenzien / Reagents used	Standard							Standard Kur	Standard Kurve / Standard curve		
Teststreifen / Antigen coated strips	SBI.DB	Ref Werte / Ref	f. Values		Gültigke	itsbereich	/ Valid	ity Range	Parameter	Α	0,032
Standardserum / Standard serum	SCI.AF	OD	0,88		OD	0,44	-	1,50		В	1,235
Negativ Kontrolle / Negative control	SCI.AE									С	4,601
Konjugat / Conjugate	SIH.CE+	Units	33,5	U/ml						D	4,135
Quantifizierungsgrenzen / Limits of quantific	U/ml	2	-	200							
Grenzwertbereich / Borderline range	U/ml	10		15							

OD Bereich / OD Range 405 nm, Standardserum / Standard serum																			
0,44		0,49	0,50 -	0,54	0,55	- 0,60	0,61	- 0,65	0,66	- 0,71	0,72	- 0,76	0,77	- 0,82	0,83	- 0,87	0,88	U/ml	Interpretation
	<	0,14	<	0,15		< 0,17		< 0,19		< 0,20		< 0,22		< 0,23		< 0,25	< 0,26	< 10,0	neg
0,14	-	0,20	0,15 -	0,23	0,17	- 0,25	0,19	- 0,28	0,20	- 0,30	0,22	- 0,33	0,23	- 0,35	0,25	- 0,38	0,26 - 0,39	10,0 - 15,0	gw / borderline
	>	0,20	>	0,23		> 0,25		> 0,28		> 0,30		> 0,33		> 0,35		> 0,38	> 0,39	> 15,0	pos

OD Bereich / OD Range 405 nm, Standardserum / Standard serum												
U/ml	0,88	0,89 - 0,96	0,97 - 1,03	1,04 - 1,11	1,12 - 1,19	1,20 - 1,27	1,28 - 1,34	1,35 - 1,42 1,43 - 1,50	Interpretation			
< 10,0	< 0,26	< 0,27	< 0,30	< 0,32	< 0,34	< 0,36	< 0,39	< 0,41 < 0,43	neg			
10,0 - 15,0	0,26 - 0,39	0,27 - 0,41	0,30 - 0,44	0,32 - 0,48	0,34 - 0,51	0,36 - 0,55	0,39 - 0,58		gw / borderline			
> 15,0	> 0,39	> 0,41	> 0,44	> 0,48	> 0,51	> 0,55	> 0,58	> 0,61 > 0,65	pos			

Formeln für spezielle Auswertesysteme Special case formulas OD =

0,447 x MV(STD) entspricht oberem cut-off/corresponds to upper cut-off

OD = **0,294** x MV(STD) entspricht unterem cut-off/corresponds to lower cut-off

 $Concentration = exp(4,601-ln(4,103/(MV(Sample)\ x0,88/\ MV(STD)-0,032)-1)/1,235)$

Institut Virion\Serion GmbH

Friedrich-Bergius-Ring 19 D-97076 Würzburg

Zusätzliche Barcodes mit Formeln für / Additional Barcodes with formulas for Revelation™ DSX / DS-Matrix ™

4PS- Formel / 4PS-formula

exp(4.601-ln(4.103/(Sample*0.880/S-0.032)-1)/1.235)



Gültigkeitsbereich / Validity Range

0.440<=S1<=1.496



If OD Sample < Parameter A

if Ti<(0.032*(S1/0.880)) then Ti=(0.032+0.001)*(S1/0.880)



If OD Sample > Parameter D

if Ti>(4.135*(S1/0.880)) then Ti=(4.135-0.001)*(S1/0.880)



If OD Negative control < Parameter A

if NC1<(0.032*(S1/0.880)) then NCi=(0.032+0.001)*(S1/0.880)



Institut Virion\Serion GmbH Friedrich-Bergius-Ring 19 D-97076 Würzburg