SERION ELISA classic

ESR1371G

CHLAMYDIA PNEUMONIAE IgG

SIH.CZ

Qualitätskontrollzertifikat / Quality Control Certificate

Kitcharge / Lot SIH.CZ IFU-Version 1371-4

Verw. Bis / Exp. 2019-09 Prüfdatum /

Date of control

10.10.2017



12

10

Verwendete Reagenzien / Reagents used	Standard							Standard Kur	Standard Kurve / Standard curve		
Teststreifen / Antigen coated strips	Ref Werte / Ref	f. Values		Gültigke	itsbereich	/ Valid	dity Range	Parameter	Α	0,024	
Standardserum / Standard serum	OD	0,84		OD	0,42	-	1,43		В	1,037	
Negativ Kontrolle / Negative control	SIH.BT									С	4,746
Konjugat / Conjugate	SHH.EA+++	Units	70,3 U	/ml						D	2,201
Quantifizierungsgrenzen / Limits of quantific	U/ml	4	-	200							
Grenzwertbereich / Borderline range	U/ml	10	-	12						•	

OD Bereich / OD Range 405 nm, Standardserum / Standard serum																			
0,42	-	0,46	0,47	- 0,52	0,53	- 0,57	0,58	- 0,62	0,63	- 0,67	0,68	- 0,73	0,74	- 0,78	0,79	- 0,83	0,84	U/ml	Interpretation
	<	0,09		< 0,11		< 0,12		< 0,13		< 0,14		< 0,15		< 0,16		< 0,17	< 0,18	< 10,0	neg
0,09	-	0,11	0,11	- 0,12	0,12	- 0,14	0,13	- 0,15	0,14	- 0,16	0,15	- 0,18	0,16	- 0,19	0,17	- 0,20	0,18 - 0,21	10,0 - 12,0	gw / borderline
	>	0,11		> 0,12		> 0,14		> 0,15		> 0,16		> 0,18		> 0,19		> 0,20	> 0,21	> 12,0	pos

OD Bereich / OD Range 405 nm, Standardserum / Standard serum												
U/ml		0,84	0,85 - 0,91	0,92 - 0,99	1,00 - 1,06	1,07 - 1,13	1,14 - 1,21	1,22 - 1,28	1,29 - 1,35 1,36 - 1,43	Interpretation		
<	10,0	< 0,18	< 0,19	< 0,20	< 0,22	< 0,24	< 0,25	< 0,27	< 0,28 < 0,30	neg		
10,0 -	12,0 12,0	, ,	0,19 - 0,22 > 0,22	0,20 - 0,24 > 0,24	0,22 - 0,26 > 0,26	0,24 - 0,28 > 0,28	0,25 - 0,29 > 0,29	0,27 - 0,31 > 0,31	0,28 - 0,33 0,30 - 0,35 > 0,33 > 0,35	gw / borderline		

Formeln für spezielle Auswertesysteme Special case formulas OD = **0,255** x MV(STD) entspricht oberem cut-off/corresponds to upper cut-off
OD = **0,219** x MV(STD) entspricht unterem cut-off/corresponds to lower cut-off

Concentration= exp(4,746-In(2,177/(MV(Sample) x0,84/ MV(STD)-0,024)-1)/1,037)

Institut Virion\Serion GmbH Friedrich-Bergius-Ring 19 D-97076 Würzburg

Zusätzliche Barcodes mit Formeln für / Additional Barcodes with formulas for Revelation™ DSX / DS-Matrix ™

4PS- Formel / 4PS-formula

exp(4.746-ln(2.177/(Sample*0.840/S-0.024)-1)/1.037)



Gültigkeitsbereich / Validity Range

0.420<=S1<=1.428



If OD Sample < Parameter A

if Ti < (0.024*(S1/0.840)) then Ti = (0.024+0.001)*(S1/0.840)



If OD Sample > Parameter D

if Ti>(2.201*(S1/0.840)) then Ti=(2.201-0.001)*(S1/0.840)



If OD Negative control < Parameter A

if NC1<(0.024*(S1/0.840)) then NCi=(0.024+0.001)*(S1/0.840)



Institut Virion\Serion GmbH Friedrich-Bergius-Ring 19 D-97076 Würzburg