SERION ELISA classic ESR106M

## **LEGIONELLA PNEUMOPHILA 1-7 IgM**

SBI.BO

## **Qualitätskontrollzertifikat / Quality Control Certificate**

Kitcharge / Lot SBI.BO IFU-Version 106-11

Verw. bis / Exp. 2019-12

07.02.2018

Prüfdatum /

**Date of control** 



Verwendete Reagenzien / Reagents used	Standard							Standard Kur	Standard Kurve / Standard curve		
Teststreifen / Antigen coated strips	SMH.DC	Ref Werte / Re	f. Values		Gültigke	itsbereich	/ Valid	lity Range	Parameter	Α	0,042
Standardserum / Standard serum	SBI.BC	10	<b>0,96</b>		OD	0,48	-	1,63		В	1,220
Negativ Kontrolle / Negative control	SBI.BB									С	5,995
Konjugat / Conjugate	SAI.CK+	Units	s <b>171 U</b> /n	nl						D	3,561
Quantifizierungsgrenzen / Limits of quantifi	U/ml	10	-	500							
Grenzwertbereich / Borderline range	U/ml	120	-	140							

OD Bereich / OD Range 405 nm, Standardserum / Standard serum																			
0,48	-	0,53	0,54	- 0,59	0,60	- 0,65	0,66	- 0,71	0,72	- 0,77	0,78	- 0,83	0,84	- 0,89	0,90 -	0,95	0,96	U/ml	Interpretation
	<	0.37		< 0.41		< 0.46		< 0.50		< 0.54		< 0.59		< 0.63	<	0.67	< 0.70	< 120,0	neg
0,37	-	0,42	0,41	- 0,47	0,46	- 0,52	0,50	- 0,57	0,54	- 0,62	0,59	- 0,67	0,63	- 0,72	0,67 -	0,77	0,70 - 0,80	120,0 - 140,0	gw / borderline
	>	0,42		> 0,47		> 0,52		> 0,57		> 0,62		> 0,67		> 0,72	>	0,77	> 0,80	> 140,0	pos

OD Bereich / OD Range 405 nm, Standardserum / Standard serum											
U/ml	0,96	0,97 - 1,04	1,05 - 1,13	1,14 - 1,21	1,22 - 1,30	1,31 - 1,38	1,39 - 1,46	1,47 - 1,55 1,56 - 1,63	Interpretation		
< 120,0	< 0,70	< 0,73	< 0,80	< 0,86	< 0,92	< 0,98	< 1,04	< 1,10 < 1,16	neg		
120,0 - 140,0 > 140,0		0,73 - 0,84 > 0,84	0,80 - 0,91 > 0,91	0,86 - 0,98 > 0,98	0,92 - 1,05 > 1,05	0,98 - 1,12 > 1,12	1,04 - 1,19 > 1.19		gw / borderline		

Formeln für spezielle Auswertesysteme Special case formulas OD = 0.838 x MV(STD) entspricht oberem cut-off/corresponds to upper cut-off OD = 0.727 x MV(STD) entspricht unterem cut-off/corresponds to lower cut-off Concentration= exp(5.995-ln(3.519/(MV(Sample) x0.96/ MV(STD)-0.042)-1)/1,22)

140 Institut Virion\Serion GmbH
120 Friedrich-Bergius-Ring 19
D-97076 Würzburg

# Zusätzliche Barcodes mit Formeln für / Additional Barcodes with formulas for Revelation™ DSX / DS-Matrix ™

#### 4PS- Formel / 4PS-formula

exp(5.995-ln(3.519/(Sample\*0.960/S-0.042)-1)/1.220)



#### Gültigkeitsbereich / Validity Range

0.480<=S1<=1.632



#### If OD Sample < Parameter A

if Ti < (0.042\*(S1/0.960)) then Ti = (0.042+0.001)\*(S1/0.960)



## If OD Sample > Parameter D

if Ti>(3.561\*(S1/0.960)) then Ti=(3.561-0.001)\*(S1/0.960)



### If OD Negative control < Parameter A

if NC1<(0.042\*(S1/0.960)) then NCi=(0.042+0.001)\*(S1/0.960)



Institut Virion\Serion GmbH Friedrich-Bergius-Ring 19 D-97076 Würzburg