**SERION ELISA classic** 

**ESR1052G** 

!New!

**HERPES SIMPLEX VIRUS 2 IgG** 

SBK.FH

## **Qualitätskontrollzertifikat / Quality Control Certificate**

Kitcharge / Lot SBK.FH IFU-Version 105-23

Verw. bis / Exp. 2020-06

07.03.2019

Prüfdatum /

Date of control



Verwendete Reagenzien / Reagents used	Standard	Standard Ku	Standard Kurve / Standard curve								
Teststreifen / Antigen coated strips	SGI.DZ	Ref Werte / Re	f. Values		Gültigke	itsbereich	/ Vali	dity Range	Parameter	Α	-0,004
Standardserum / Standard serum	SKI.BB	O	0,91		OD	0,46	-	1,55		В	0,936
Negativ Kontrolle / Negative control	SKI.BA									С	5,814
Konjugat / Conjugate	SMI.CU+++	Units	s <b>220 U/m</b>	I						D	2,267
Quantifizierungsgrenzen / Limits of quantifi	U/ml	10	-	1000							
Grenzwertbereich / Borderline range	U/ml	20	-	30							

OD Bereich / OD Range 405 nm, Standardserum / Standard serum																			
0,46	-	0,50	0,51	- 0,56	0,57	- 0,62	0,63	- 0,67	0,68	- 0,73	0,74	- 0,79	0,80	- 0,84	0,85	- 0,90	0,91	U/ml	Interpretation
	<	0,08		< 0,09		< 0,10		< 0,11		< 0,12		< 0,13		< 0,14		< 0,14	< 0,15	< 20	),0 neg
0,08	-	0,11	0,09	- 0,12	0,10	- 0,14	0,11	- 0,15	0,12	- 0,16	0,13	- 0,18	0,14	- 0,19	0,14	- 0,20	0,15 - 0,21	20,0 - 30	),0 gw / borderline
	>	0,11		> 0,12		> 0,14		> 0,15		> 0,16		> 0,18		> 0,19		> 0,20	> 0,21	> 30	),0 pos

OD Bereich / OD Range 405 nm, Standardserum / Standard serum													
U/ml			0,91	0,92 - 0,99	1,00 - 1,07	1,08 - 1,15	1,16 - 1,23	1,24 - 1,31	1,32 - 1,39	1,40 - 1,47 1,48 - 1,55	Interpretation		
	<	20,0	< 0,15	< 0,16	< 0,17	< 0,18	< 0,20	< 0,21	< 0,22	< 0,24 < 0,25	neg		
20,0	- >	30,0 30,0	0,15 - 0,21 > 0,21	0,16 - 0,22 > 0,22	0,17 - 0,24 > 0,24	0,18 - 0,26 > 0,26	0,20 - 0,28 > 0,28	0,21 - 0,29 > 0,29	0,22 - 0,31 > 0,31		gw / borderline pos		

Formeln für spezielle Auswertesysteme Special case formulas OD = **0,232** x MV(STD) entspricht oberem cut-off/corresponds to upper cut-off
OD = **0,162** x MV(STD) entspricht unterem cut-off/corresponds to lower cut-off

Concentration= exp(5,814-In(2,271/(MV(Sample) x0,91/ MV(STD)+0,004)-1)/0,936)

30 Institut Virion\Serion GmbH20 Friedrich-Bergius-Ring 19

Friedrich-Bergius-Ring 19 D-97076 Würzburg

# Zusätzliche Barcodes mit Formeln für / Additional Barcodes with formulas for Revelation™ DSX / DS-Matrix ™

#### 4PS- Formel / 4PS-formula

exp(5.814-ln(2.271/(Sample\*0.910/S+0.004)-1)/0.936)



#### Gültigkeitsbereich / Validity Range

0.455<=S1<=1.547



#### If OD Sample < Parameter A

if Ti<(-0.004\*(S1/0.910)) then Ti=(-0.004+0.001)\*(S1/0.910)



# If OD Sample > Parameter D

if Ti>(2.267\*(S1/0.910)) then Ti=(2.267-0.001)\*(S1/0.910)



### If OD Negative control < Parameter A

if NC1<(-0.004\*(S1/0.910)) then NCi=(-0.004+0.001)\*(S1/0.910)



Institut Virion\Serion GmbH Friedrich-Bergius-Ring 19 D-97076 Würzburg