SERION ELISA classic ESR118M

HELICOBACTER PYLORI IgM

SFI.EM

Qualitätskontrollzertifikat / Quality Control Certificate

Kitcharge / Lot SFI.EM IFU-Version 118-14

Verw. bis / Exp. 2020-05

26.06.2018

Prüfdatum /

Date of control



30

20

Verwendete Reagenzien / Reagents used	Standard							Standard Kur	Standard Kurve / Standard curve			
Teststreifen / Antigen coated strips	SEI.BI	Ref Werte / Ref	. Value	es		Gültigkei	itsbereich	/ Vali	dity Range	Parameter	Α	0,020
Standardserum / Standard serum	SFI.BK	OD	0,82			OD	0,41	-	1,39		В	1,026
Negativ Kontrolle / Negative control	SFI.BI										С	5,078
Konjugat / Conjugate	SEI.BR+	Units	46,8	U/ml							D	3,655
Quantifizierungsgrenzen / Limits of quantifi	U/ml	10) .		300							
Grenzwertbereich / Borderline range	U/ml	20) .		30	<u> </u>						

OD Bereich / OD Range 405 nm, Standardserum / Standard serum																			
0,41	-	0,45	0,46 -	0,50	0,51	- 0,55	0,56	- 0,61	0,62	- 0,66	0,67	- 0,71	0,72	- 0,76	0,77	- 0,81	0,82	U/ml	Interpretation
	<	0,21	<	0,24		< 0,26		< 0,29		< 0,31		< 0,34		< 0,36		< 0,39	< 0,40	< 20,0	neg
0,21	-	0,30	0,24 -	0,33	0,26	- 0,37	0,29	- 0,41	0,31	- 0,44	0,34	- 0,48	0,36	- 0,51	0,39	- 0,55	0,40 - 0,57	20,0 - 30,0	gw / borderline
	>	0,30	>	- 0,33		> 0,37		> 0,41		> 0,44		> 0,48		> 0,51		> 0,55	> 0,57	> 30,0	pos

OD Bereich / OD Range 405 nm, Standardserum / Standard serum											
U/ml	0,82	0,83 - 0,89	0,90 - 0,96	0,97 - 1,04	1,05 - 1,11	1,12 - 1,18	1,19 - 1,25	1,26 - 1,32 1,33 - 1,39	Interpretation		
< 20,0	< 0,40	< 0,42	< 0,45	< 0,49	< 0,52	< 0,56	< 0,59	< 0,63 < 0,66	neg		
20,0 - 30,0 > 30,0	0,40 - 0,57 > 0,57	0,42 - 0,60 > 0,60	0,45 - 0,65 > 0,65	0,49 - 0,70 > 0,70	0,52 - 0,75 > 0,75	0,56 - 0,80 > 0,80	0,59 - 0,85 > 0,85		gw / borderline pos		

Formeln für spezielle Auswertesysteme Special case formulas $\begin{array}{lll} \text{OD} = & \textbf{0,697} & \text{x MV(STD)} \text{ entspricht oberem cut-off/corresponds to upper cut-off} \\ \text{OD} = & \textbf{0,493} & \text{x MV(STD)} \text{ entspricht unterem cut-off/corresponds to lower cut-off} \\ \end{array}$

 $Concentration = \exp(5,078 - \ln(3,635/(\text{MV(Sample}) \ x0,82/\ \text{MV(STD)} - 0,02) - 1)/1,026)$

Institut Virion∖Serion GmbH Friedrich-Bergius-Ring 19 D-97076 Würzburg

Zusätzliche Barcodes mit Formeln für / Additional Barcodes with formulas for Revelation™ DSX / DS-Matrix ™

4PS- Formel / 4PS-formula

exp(5.078-ln(3.635/(Sample*0.820/S-0.020)-1)/1.026)



Gültigkeitsbereich / Validity Range

0.410<=S1<=1.394



If OD Sample < Parameter A

if Ti<(0.020*(S1/0.820)) then Ti=(0.020+0.001)*(S1/0.820)



If OD Sample > Parameter D

if Ti>(3.655*(S1/0.820)) then Ti=(3.655-0.001)*(S1/0.820)



If OD Negative control < Parameter A

if NC1<(0.020*(S1/0.820)) then NCi=(0.020+0.001)*(S1/0.820)



Institut Virion\Serion GmbH Friedrich-Bergius-Ring 19 D-97076 Würzburg