SERION ELISA classic

ESR1052G

HERPES SIMPLEX VIRUS 2 IgG

SEI.CC

Qualitätskontrollzertifikat / Quality Control Certificate

Kitcharge / Lot SEI.CC IFU-Version 105-22

Verw. bis / Exp. 2020-04

24.05.2018

Prüfdatum /

Date of control



30

20

Verwendete Reagenzien / Reagents used	Standard							Standard Kui	Standard Kurve / Standard curve		
Teststreifen / Antigen coated strips	SCI.CY	Ref Werte / Ref	f. Values		Gültigke	itsbereich	/ Valid	lity Range	Parameter	Α	-0,014
Standardserum / Standard serum	SDI.DH	OD	0,88		OD	0,44	-	1,50		В	0,888
Negativ Kontrolle / Negative control	SDI.DG									С	5,793
Konjugat / Conjugate	SDI.EB+++	Units	207 U/n	nl						D	2,227
Quantifizierungsgrenzen / Limits of quantific	U/ml	10	-	1000							
Grenzwertbereich / Borderline range	U/ml	20	-	30							

OD Bereich / OD Range 405 nm, Standardserum / Standard serum																			
0,44	-	0,49	0,50	- 0,54	0,55	- 0,60	0,61	- 0,65	0,66	- 0,71	0,72	- 0,76	0,77	- 0,82	0,83	- 0,87	0,88	U/ml	Interpretation
	<	0.08		< 0.09	1	< 0.10		< 0.11		< 0.12		< 0.13		< 0.14		< 0.15	< 0.16	< 20.0	neg
0,08	-		0,09	- 0,14	0,10	- 0,15	0,11	- 0,16	0,12	- 0,18	0,13	- 0,19	0,14	- 0,21	0,15	- 0,22	0,16 - 0,23		gw / borderline
	>	0,12		> 0,14		> 0,15		> 0,16		> 0,18		> 0,19		> 0,21		> 0,22	> 0,23	> 30,0	pos

OD Bereich / OD Range 405 nm, Standardserum / Standard serum												
U/ml		0,88	0,89 - 0,96	0,97 - 1,03	1,04 - 1,11	1,12 - 1,19	1,20 - 1,27	1,28 - 1,34	1,35 - 1,42 1,43 - 1,50	Interpretation		
	< 20,0	< 0,16	< 0,17	< 0,18	< 0,20	< 0,21	< 0,22	< 0,24	< 0,25 < 0,27	neg		
20,0	- 30,0 > 30,0	0,16 - 0,23 > 0,23	0,17 - 0,24 > 0,24	0,18 - 0,26 > 0,26	0,20 - 0,28 > 0,28	0,21 - 0,30 > 0,30	0,22 - 0,32 > 0,32	0,24 - 0,34 > 0,34	0,25 - 0,36 0,27 - 0,38 > 0,36 > 0,38	gw / borderline pos		

Formeln für spezielle Auswertesysteme Special case formulas $\begin{array}{lll} \text{OD} = & \textbf{0,256} & \text{x MV(STD)} \text{ entspricht oberem cut-off/corresponds to upper cut-off} \\ \text{OD} = & \textbf{0,180} & \text{x MV(STD)} \text{ entspricht unterem cut-off/corresponds to lower cut-off} \\ \end{array}$

 $Concentration = exp(5,793-In(2,241/(MV(Sample)\ x0,88/\ MV(STD)+0,014)-1)/0,888)$

Institut Virion∖Serion GmbH Friedrich-Bergius-Ring 19 D-97076 Würzburg

Zusätzliche Barcodes mit Formeln für / Additional Barcodes with formulas for Revelation™ DSX / DS-Matrix ™

4PS- Formel / 4PS-formula

exp(5.793-ln(2.241/(Sample*0.880/S+0.014)-1)/0.888)



Gültigkeitsbereich / Validity Range

0.440<=S1<=1.496



If OD Sample < Parameter A

if Ti < (-0.014*(S1/0.880)) then Ti = (-0.014+0.001)*(S1/0.880)



If OD Sample > Parameter D

if Ti>(2.227*(S1/0.880)) then Ti=(2.227-0.001)*(S1/0.880)



If OD Negative control < Parameter A

if NC1<(-0.014*(S1/0.880)) then NCi=(-0.014+0.001)*(S1/0.880)



Institut Virion\Serion GmbH Friedrich-Bergius-Ring 19 D-97076 Würzburg