# **SERION ELISA classic**

**ESR1372M** 

## **CHLAMYDIA TRACHOMATIS IgM**

# SAI.BR

#### **Qualitätskontrollzertifikat / Quality Control Certificate**

Kitcharge / Lot SAI.BR IFU-Version 1372-6

12.01.2018 Prüfdatum /

Verw. bis / Exp. 2019-09

Date of control



14

10

Verwendete Reagenzien / Reagents used	Standard		Standard Kur	Standard Kurve / Standard curve							
Teststreifen / Antigen coated strips	Ref Werte / Ref. Values			Gültigkeitsbereich / Validity Range				Parameter	Α	0,010	
Standardserum / Standard serum	SIH.BS	OD	0,93		OD	0,47	-	1,58		В	0,885
Negativ Kontrolle / Negative control	SIH.BR									С	5,288
Konjugat / Conjugate	SKH.AQ++	Units	31,1 L	J/ml						D	5,661
Quantifizierungsgrenzen / Limits of quantific	U/ml	5	-	180							
Grenzwertbereich / Borderline range	U/ml	10		14							

OD Bereich / OD Range 405 nm, Standardserum / Standard serum																				
0,47	-	0,51	0,52	- 0,57	0,58	- 0,63	0,64	- 0,69	0,70	- 0,75	0,76	- 0,80	0,81	- 0,86	0,87	- 0,92	0,93	U/mI		Interpretation
	<	0,21		< 0,23		< 0,25		< 0,28		< 0,30		< 0,33		< 0,35		< 0,38	< 0,39	<	10,0	neg
0,21	-	0,26	0,23	- 0,29	0,25	- 0,33	0,28	- 0,36	0,30	- 0,39	0,33	- 0,42	0,35	- 0,45	0,38	- 0,48	0,39 - 0,50	10,0 -	14,0	gw / borderline
	>	0,26		> 0,29		> 0,33		> 0,36		> 0,39		> 0,42		> 0,45		> 0,48	> 0,50	>	14,0	pos

OD Bereich / OD Range 405 nm, Standardserum / Standard serum												
U/ml		0,93	0,94 - 1,01	1,02 - 1,09	1,10 - 1,17	1,18 - 1,26	1,27 - 1,34	1,35 - 1,42	1,43 - 1,50 1,51 - 1,	58 Interpretation		
<	10,0	< 0,39	< 0,41	< 0,44	< 0,48	< 0,51	< 0,55	< 0,58	< 0,61 < 0,	65 neg		
10,0 -	14,0	, ,	0,41 - 0,52	0,44 - 0,57	0,48 - 0,61	, ,	0,55 - 0,70	0,58 - 0,74		83 gw / borderline		
>	14,0	> 0,50	> 0,52	> 0,57	> 0,61	> 0,66	> 0,70	> 0,74	> 0,79 > 0,	83 pos		

Formeln für spezielle Auswertesysteme Special case formulas OD = **0,543** x MV(STD) entspricht oberem cut-off/corresponds to upper cut-off
OD = **0,415** x MV(STD) entspricht unterem cut-off/corresponds to lower cut-off

Concentration= exp(5,288-In(5,651/(MV(Sample) x0,93/ MV(STD)-0,01)-1)/0,885)

Institut Virion\Serion GmbH Friedrich-Bergius-Ring 19 D-97076 Würzburg

# Zusätzliche Barcodes mit Formeln für / Additional Barcodes with formulas for Revelation™ DSX / DS-Matrix ™

#### 4PS- Formel / 4PS-formula

exp(5.288-ln(5.651/(Sample\*0.930/S-0.010)-1)/0.885)



#### Gültigkeitsbereich / Validity Range

0.465<=S1<=1.581



#### If OD Sample < Parameter A

if Ti < (0.010\*(S1/0.930)) then Ti = (0.010+0.001)\*(S1/0.930)



## If OD Sample > Parameter D

if Ti>(5.661\*(S1/0.930)) then Ti=(5.661-0.001)\*(S1/0.930)



### If OD Negative control < Parameter A

if NC1<(0.010\*(S1/0.930)) then NCi=(0.010+0.001)\*(S1/0.930)



Institut Virion\Serion GmbH Friedrich-Bergius-Ring 19 D-97076 Würzburg