## SERION ELISA classic

ESR139M

## **CAMPYLOBACTER JEJUNI IgM**

SAI.DI

### **Qualitätskontrollzertifikat / Quality Control Certificate**

Kitcharge / Lot SAI.DI IFU-Version 139-10

17.01.2018 **Prüfdatum** /

Verw. bis / Exp. 2019-11

Date of control



60

40

Verwendete Reagenzien / Reagents used	Standard							Standard Kur	Standard Kurve / Standard curve		
Teststreifen / Antigen coated strips	SLH.DN	Ref Werte / Re	f. Values		Gültigke	itsbereich	/ Valid	lity Range	Parameter	Α	-0,002
Standardserum / Standard serum	SAI.BI	00	0,86		OD	0,43	-	1,46		В	0,825
Negativ Kontrolle / Negative control	SAI.BH									С	7,167
Konjugat / Conjugate	SFH.AB++	Units	s <b>348 U</b> /r	nl						D	3,413
Quantifizierungsgrenzen / Limits of quantific	U/ml	20	-	2000							
Grenzwertbereich / Borderline range	U/ml	40	-	60							

OD Bereich / OD Range 405 nm, Standardserum / Standard serum																			
0,43	-	0,47	0,48	- 0,53	0,54	- 0,58	0,59	- 0,64	0,65	- 0,69	0,70	- 0,74	0,75	- 0,80	0,81	- 0,85	0,86	U/ml	Interpretation
	<	0,09		< 0,11		< 0,12		< 0,13		< 0,14		< 0,15		< 0,16		< 0,17	< 0,18	< 40,0	neg
0,09	-	0,13	0,11	- 0,15	0,12	- 0,16	0,13	- 0,18	0,14	- 0,19	0,15	- 0,21	0,16	- 0,23	0,17	- 0,24	0,18 - 0,25	40,0 - 60,0	gw / borderline
	>	0,13		> 0,15		> 0,16		> 0,18		> 0,19		> 0,21		> 0,23		> 0,24	> 0,25	> 60,0	pos

OD Bereich / OD Range 405 nm, Standardserum / Standard serum											
U/mI	0,86	0,87 - 0,94	0,95 - 1,01	1,02 - 1,09	1,10 - 1,16	1,17 - 1,24	1,25 - 1,31	1,32 - 1,39 1,40 - 1,46	Interpretation		
< 40,0	< 0,18	< 0,19	< 0,20	< 0,22	< 0,24	< 0,25	< 0,27	< 0,28 < 0,30	neg		
40,0 - 60,0 > 60,0	0,18 - 0,25 > 0,25	0,19 - 0,26 > 0,26	0,20 - 0,28 > 0,28	0,22 - 0,31 > 0,31	0,24 - 0,33 > 0,33	0,25 - 0,35 > 0,35	0,27 - 0,37 > 0,37		gw / borderline		

Formeln für spezielle Auswertesysteme Special case formulas OD = **0,289** x MV(STD) entspricht oberem cut-off/corresponds to upper cut-off

 $\label{eq:ode_ode_ode_ode_ode} OD = \qquad \textbf{0,211} \qquad x \ \text{MV(STD)} \ \text{entspricht unterem cut-off/corresponds to lower cut-off}$ 

 $Concentration = exp(7,167-In(3,415/(MV(Sample)\ x0,86/\ MV(STD)+0,002)-1)/0,825)$ 

Institut Virion\Serion GmbH Friedrich-Bergius-Ring 19 D-97076 Würzburg

# Zusätzliche Barcodes mit Formeln für / Additional Barcodes with formulas for Revelation™ DSX / DS-Matrix ™

#### 4PS- Formel / 4PS-formula

exp(7.167-ln(3.415/(Sample\*0.860/S+0.002)-1)/0.825)



#### Gültigkeitsbereich / Validity Range

0.430<=S1<=1.462



#### If OD Sample < Parameter A

if Ti < (-0.002\*(S1/0.860)) then Ti = (-0.002+0.001)\*(S1/0.860)



## If OD Sample > Parameter D

if Ti>(3.413\*(S1/0.860)) then Ti=(3.413-0.001)\*(S1/0.860)



## If OD Negative control < Parameter A

if NC1 < (-0.002\*(S1/0.860)) then NCi = (-0.002+0.001)\*(S1/0.860)



Institut Virion\Serion GmbH Friedrich-Bergius-Ring 19 D-97076 Würzburg