SERION ELISA classic

ESR1052M

HERPES SIMPLEX VIRUS 2 IgM

SDI.BT

Qualitätskontrollzertifikat / Quality Control Certificate

Kitcharge / Lot SDI.BT IFU-Version 105-22

Prüfdatum /

Verw. bis / Exp. 2020-04

i raidataiii /

13.04.2018

Date of control



30

20

Verwendete Reagenzien / Reagents used	Standard							Standard Kurve / Standard curve			
Teststreifen / Antigen coated strips	SAI.FL	Ref Werte / Re	f. Values		Gültigke	itsbereich	/ Valid	dity Range	Parameter	Α	0,107
Standardserum / Standard serum	SDI.AO	00	0,77		OD	0,39	-	1,31		В	1,163
Negativ Kontrolle / Negative control	SDI.AN									С	3,536
Konjugat / Conjugate	SBI.EN++	Units	s 10,4 U /	ml						D	3,432
Quantifizierungsgrenzen / Limits of quantific	U/mI	10	-	500							
Grenzwertbereich / Borderline range	U/ml	20	-	30							

OD Bereich / OD Range 405 nm, Standardserum / Standard serum																		
0,39	-	0,42	0,43	- 0,47	0,48	- 0,52	0,53	- 0,57	0,58	- 0,62	0,63 -	0,66	0,67 - 0,71	0,72	- 0,76	0,77	U/ml	Interpretation
	<	0,66		< 0,74		< 0,82		< 0,90		< 0,98	<	1,06	< 1,14		< 1,22	< 1,26	< 20,0	neg
0,66	-	0,86	0,74	- 0,96	0,82	- 1,06	0,90	- 1,17	0,98	- 1,27	1,06 -	1,37	1,14 - 1,47	1,22	- 1,58	1,26 - 1,64	20,0 - 30,0	gw / borderline
	>	0,86		> 0,96		> 1,06		> 1,17		> 1,27	>	1,37	> 1,47		> 1,58	> 1,64	> 30,0	pos

OD Bereich / OD Range 405 nm, Standardserum / Standard serum											
U/ml	0,77	0,78 - 0,84	0,85 - 0,90	0,91 - 0,97	0,98 - 1,04	1,05 - 1,11	1,12 - 1,17	1,18 - 1,24 1,25 - 1,31	Interpretation		
< 20,0	< 1,26	< 1,33	< 1,44	< 1,55	< 1,66	< 1,77	< 1,88	< 1,99 < 2,10	neg		
20,0 - 30,0 > 30.0	1,26 - 1,64 > 1,64	1,33 - 1,72 > 1,72	1,44 - 1,86 > 1,86	1,55 - 2,01 > 2,01	1,66 - 2,15 > 2,15	1,77 - 2,29 > 2.29	1,88 - 2,44 > 2,44	1,99 - 2,58 2,10 - 2,73 > 2,58 > 2,73	gw / borderline		

Formeln für spezielle Auswertesysteme Special case formulas

 $Concentration = exp(3,536-ln(3,325/(MV(Sample)\ x0,77/\ MV(STD)-0,107)-1)/1,163)$

Institut Virion\Serion GmbH Friedrich-Bergius-Ring 19 D-97076 Würzburg

Zusätzliche Barcodes mit Formeln für / Additional Barcodes with formulas for Revelation™ DSX / DS-Matrix ™

4PS- Formel / 4PS-formula

exp(3.536-ln(3.325/(Sample*0.770/S-0.107)-1)/1.163)



Gültigkeitsbereich / Validity Range

0.385<=S1<=1.309



If OD Sample < Parameter A

if Ti<(0.107*(S1/0.770)) then Ti=(0.107+0.001)*(S1/0.770)



If OD Sample > Parameter D

if Ti>(3.432*(S1/0.770)) then Ti=(3.432-0.001)*(S1/0.770)



If OD Negative control < Parameter A

if NC1<(0.107*(S1/0.770)) then NCi=(0.107+0.001)*(S1/0.770)



Institut Virion\Serion GmbH Friedrich-Bergius-Ring 19 D-97076 Würzburg